



L I G H T & V I S I O N




VISION

- LARGER FIELD OF VIEW
- SUPERIOR DEPTH FOCUS
- CLEAR & BRIGHT VISION OVER THE FULL FIELD

COMFORT

- LIGHTER & SMALLER TELESCOPES
- NO DISTORTION & FATIGUE
- ANTI-ALLERGIC **TITANIUM** FRAME
- ALUMINIUM** FRAME = LIGHTER **NEW**
- ERGONOMICAL
- New@ForYouMed



CONTACT INFO:

ForYouMed
 medical services & supplies

www.ForYouMed.com


LIGHT



TOTAL WEIGHT
100 gr

FLEXIBLE HEAD UNIT - INCL SOFT HEADBAND

LARGER AND UNIFORM




UNIFORM BRIGHTNESS

NO SHADOWS

Brightness (lux) & Operating Time (hr)

Diameter of illumination @WD400mm		
LLM1582	LLM1574	LLM1566
Ø82 mm	Ø74 mm	Ø66 mm

Weight op LED Lens Module		
LLM1582	LLM1574	LLM1566
35.1 g	32.8 g	34.2 g



Dr. Bernard Depypere

Prof. Dr. Koen Van Landuyt

TESTIMONIAL

"The brightness, longevity and lightweight of the loupes and ledlights are amazing! With an attractive price, on site support and a great service."

Plastic surgeon
Prof. Dr. Alexis Verpaele
Ghent, Belgium



EXCLUSIVE
DISTRIBUTOR



VISION

"AS A SURGEON OR PROFESSIONAL
NOTHING BUT THE BEST IS GOOD ENOUGH,
FOR YOU..."

WE OFFER THE ULTIMATE VIEWING EXPERIENCE DUE TO OUR CUSTOM MADE LOUPES WHICH ARE DEVELOPED THROUGH CONTINUOUS R&D TO ACHIEVE PERFECTION IN HIGH OPTICAL TECHNOLOGY.

TITANIUM



TTL Loupes comparison table*

Magnification:	2.5x	2.8x	3.2x	3.5x	5x
Field of vision:	160 mm	130 mm	86 mm	98 mm	62 mm
Weight:	32.8 g	40 g	56.6 g	73.9 g	75.4 g

ALUMINIUM

20 gr



TTL Loupes comparison table*

Magnification:	2.5x	2.8x	3.2x	3.5x	5x
Field of vision:	160 mm	130 mm	86 mm	98 mm	62 mm
Weight:	20 g	27.2 g	43.8 g	61.1 g	62.6 g

* Field of vision is given for an example 40 cm working distance.



PRODUCT : LOOKS + MAG 3.5X

NO DISTORTION

LIGHT WEIGHT

COMPACT

EXTENDED FIELD