



L I G H T & V I S I O N




VISION

- LARGER FIELD OF VIEW
- SUPERIOR DEPTH FOCUS
- CLEAR & BRIGHT VISION OVER THE FULL FIELD

COMFORT

- LIGHTER & SMALLER TELESCOPES
- NO DISTORTION & FATIGUE
- ANTI-ALLERGIC **TITANIUM** FRAME
- ALUMINIUM** FRAME = LIGHTER **NEW**
- ERGONOMICAL
- New@ForYouMed



CONTACT INFO:

ForYouMed
 medical services & supplies

www.ForYouMed.com

LIGHT




CLARITY

UNIFORM BRIGHTNESS

NO SHADOWS


LARGER AND UNIFORM



Brightness (lux) & Operating Time (hr)

Diameter of illumination @WD400mm		
LLM1582	LLM1574	LLM1566
Ø82 mm	Ø74 mm	Ø66 mm

Weight op LED Lens Module		
LLM1582	LLM1574	LLM1566
35.1 g	32.8 g	34.2 g



TESTIMONIAL

" The brightness, longevity and lightweight of the loupes and ledlights are amazing! With an attractive price, on site support and a great service "

Dentist **Pieter Jan Vandecasteele, Ghent**

" In my opinion loupes are an essential part of a dentists' equipment if you want to give your patients high-end state of the art dental treatment. Therefore I would highly recommend these loupes. They are lightweight and comfortable to wear, and the service and support are excellent "

Dentist **Gilles D'Haese, Ghent**



EXCLUSIVE DISTRIBUTOR



VISION

"AS A DENTIST, HYGIENIST OR ORTHODONTIST, NOTHING BUT THE BEST IS GOOD ENOUGH, FOR YOU..."

WE OFFER THE ULTIMATE VIEWING EXPERIENCE DUE TO OUR CUSTOM MADE LOUPES WHICH ARE DEVELOPED THROUGH CONTINUOUS R&D TO ACHIEVE PERFECTION IN HIGH OPTICAL TECHNOLOGY.

TITANIUM



TTL Loupes comparison table*

Magnification:	2.5x	2.8x	3.2x	3.5x	5x
Field of vision:	160 mm	130 mm	86 mm	98 mm	62 mm
Weight:	32.8 g	40 g	56.6 g	73.9 g	75.4 g

ALUMINIUM

20 gr



TTL Loupes comparison table*

Magnification:	2.5x	2.8x	3.2x	3.5x	5x
Field of vision:	160 mm	130 mm	86 mm	98 mm	62 mm
Weight:	20 g	27.2 g	43.8 g	61.1 g	62.6 g

* Field of vision is given for an example 40 cm working distance.

NO DISTORTION

LIGHT WEIGHT

COMPACT

EXTENDED FIELD



PRODUCT : LOOKS + MAG 3.5X